

## 2024/7/24 (Wed)

13:00-20:30 (JST)

### Experiment Supervisor:

Koichi Yano

### Participants:

- Shiori Kajikawa (all)
- Shoya Inoue (all)
- Kei Hato (-15:00)

### Experiment:

- Lab safety guide
- gDNA purification of *B. subtilis* natto BEST195
- planting 8 different strains of *B. subtilis*, *B. subtilis* subsp. Natto, *E. Coli*

### Results:

- Participants learned the safety rules in Sue'tsugu lab.
- gDNA precipitation was visually confirmed

### Additional Notes:

- Participants must make sure to take note of, and share cautions within the lab with other members

## 2024/07/25 (Thu)

10:00-18:30 (JST)

### Experiment Supervisor:

Koichi Yano

### Participants:

- Mizuho Sakai (11:30-)
- Rikuto Fukushima (-16:30)
- Shiori Kajikawa (-13:00)
- Kei Hato (11:30-15:00)

### Experiment

- PCR of *aprE* fragment
- TAE-EtBr agarose gel creation (2% and 0.8%)

- PCR of extracted gDNA
- Moving bacteria to freezing glycerol stock

## **Results**

None

# 2024/07/26 (Fri)

12:30-18:30 (JST)

## **Experiment Supervisor:**

Koichi Yano

## **Participants:**


- Mizuho Sakai (all)
- Rikuto Fukushima(13:00-)
- Shiori Kajikawa (15:50-)

## **Experiment:**

- Autoclaving effluent
- Electrophoresis yesterday's PCR product
- PCR of DNA extracted yesterday through different cycle setting
- Researching about BioBrick RFC10

## **Results:**

- No bands except lane 10 (plasmid) or lane 1/12(ladder)showed.
- Upon troubleshooting, we found that the PCR setup, conducted yesterday, was done in 2-steps instead of 3



#### Additional Notes:

2024/07/27 (Sat)

10:30-18:00 (JST)

#### Experiment Supervisor:

Kazuyuki Fujimitsu


#### Participants:

- Rikuto Fukushima (-15:30)
- Shiori Kajikawa (11:30-)
- Tatsuhiko Akiyama (17:00-)
- Yukiya Horiba (17:00-)


#### Experiment:

- Electrophoresis of yesterday's PCR product
- Troubleshooting
- Extraction of gDNA from *B. subtilis* natto BEST195 (This time, use Wizard kit)
- PCR and electrophoresis (3rd time) of PCR product
- Pipette training session for Tatsuhiko and Yukiya

#### Results:



- No bands except ladders were seen in electrophoresis.
- Made alterations to PCR= increase extension time, lower annealing temperature, increase the addition of primers and DNA sample, reconducted DNA extraction which may have failed first time



# 2024/07/29 (Mon)

10:30- 19:30 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Rikuto Fukushima (all)
- Shiori Kajikawa (-14:30)
- Kei Hato (-17:00)
- Saki Tsuchiya (13:30-18:00)
- Mizuho Sakai (14:30-15:30)

**Experiment:**

<Wet lab>

- Clarification of plans for this week (increased production of NK and fibrin plate assay)
- Ordering necessary material for mainly this week's experiments (fibrin, thrombin)
- PCR and electrophoresis of vector plasmid


<Education>

- PCR trial for lab session on Aug 4th

**Results:**

<Wet lab>

- Experiment was planned for this week
- PCR confirmed via electrophoresis



#### <Education>

•PCR & electrophoresis committed 3 times, all negative results. Resume trial tomorrow, but adding water and heating in microwave when breaking down the soybeans

## 2024/07/30 (Tue)

12:30-17:00 (JST)

#### Experiment Supervisor:

Koichi Yano

#### Participants:

- Shiori Kajikawa (all)
- Saki Tsuchiya (13:00-17:00)
- Mizuho Sakai (13:00-)
- Kei Hato (13:30-)
- Shoya Inoue (14:00-)

#### Experiment:

##### <Wet lab>

- Creation of culture medium for *B. subtilis*, *B. subtilis* Natto, and preculture of *E. Coli* C600
- Ordering primers for colony PCR fragments

<Education>

PCR trial for experiment session on Aug. 4

**Results:**

<Wet lab>

None

<Education>

•PCR failed. No bands except the ladder.

## 2024/07/31 (Wed)

11:00-18:00 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Shiori Kajikawa (all)
- Mizuho Sakai (11:30-)
- Rikuto Fukushima (14:00-)

**Experiment:**

- Preparation of E. coli C600 competent cell
- Creation of LB medium for transformation

**Results:**

- At the end, we realized that the “E. coli” cells we were using was actually B. subtilis NEST115.
- Thus, we made a glycerol stock of competent B. subtilis NEST115 cells.
- We will re-conduct the same experiment using E. coli C600

**Additional Notes:**

## 2024/08/01 (Thu)

10:00-17:00 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Shiori Kajikawa (all)
- Kei Hato (all)

- Rikuto Fukushima (13:00-)
- Ryuzo Kijima (16:00-)

**Experiment:**

- Preparation of E. Coli C600 competent cell glycerol stock

**Results:**

- Made a glycerol stock of competent E. Coli C600 competent cell

**Additional Notes:**

- Creation of phosphate buffer

## 2024/08/02 (Fri)

11:00-15:00 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Shiori Kajikawa (all)
- Shoya Inoue (-13:00)
- Rikuto Fukushima (13:30-)
- Ryuzo Kijima (14:00-)

**Experiment:**

- DNA Column purification of aprE and plasmid fragment, DNA quantity determination using Nanodrop and Quantus

**Results:**

- NanoDrop DNA concentration measurement was 168.8ng/μl, 141.8ng/μl, 137.8 ng/μl, and 103.5 ng/μl for tubes A, B, C, D respectively, for the PCR samples from 7/27 and 7/29
- Quantus DNA concentration measurement was 139.7nM, 140.5nM, 134.8nM, 22.6nM for tubes A, B, C, D respectively

## 2024/08/05 (Mon)

10:00-21:00 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**



- Kei Hato (-13:00)
- Misaki Ozawa (11:00-18:30)
- Ryuzo Kijima (13:30-18:30)
- Shiori Kajikawa (13:30-)

#### Experiment:

- infusion cloning of aprE fragments with pNK1
- transformation into E. coli DH5α
- Competence measurement of E. coli C600 with transformation

#### Results:

None

## 2024/08/06 (Tue)

11:00-19:30 (JST)

#### Experiment Supervisor:

Koichi Yano

#### Participants:

- Shiori Kajikawa (-14:00)
- Yukiya Horiba (13:30-)

#### Experiment:



- Colony PCR of transformed E. Coli DH5α/pNK1 and electrophoresis




Tube No.	Template	
1-23	DH5α infusion colonies	
24	Plasmid <05> Yano	
25	DH5α/<05> colony	
PCR program		
predenaturation	98°C, 3 min	
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	130 cycles
extension	68°C, 60sec	
final extension	68°C, 7 min	
hold	12°C	

- Inoculation of successful colonies to liquid medium

#### Results:

- Colonies were observed on the transformation plate
- All colonies were estimated to be about the targeted length (7.5kbp)





## 2024/08/07 (Wed)

10:00-17:30 (JST)

### **Experiment Supervisor:**

Koichi Yano

### **Participants:**

- Kei Hato (all)
- Misaki Ozawa (11:00-17:00)
- Rikuto Fukushima (13:00-15:00)

### **Experiment:**

- Trial of creating milk plate/fibrin plate assay

### **Results:**

None

## 2024/08/08 (Thu)

10:00-20:00 (JST)

### **Experiment Supervisor:**

Koichi Yano

### **Participants:**

- Kei Hato (10:00-17:30)
- Rikuto Fukushima
- Misaki Ozawa (14:00-20:00)

### **Experiment:**


- Creation of milk plate/fibrin plate assay
- PCR and electrophoresis for miniprep pNK1(DH5α)

Tube No.	Template DNA	NanoDrop (ng/ $\mu$ L)
1	pNK1(DH5 $\alpha$ ) #4 240807	
2	pNK1(DH5 $\alpha$ ) #8 240807	
3	pNK1(DH5 $\alpha$ ) #14 240807	
4	pNK1(DH5 $\alpha$ ) #19 240807	

PCR program		
predenaturation	98°C, 3 min	
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	730 cycles
extension	68°C, 30sec	
final extension	68°C, 7 min	
hold	12°C	

## Results:



## Additional Notes:

2024/08/09 (Fri)

10:00-18:00 (JST)

## Experiment Supervisor:

Koichi Yano

## Participants:



• Shiori Kajikawa (all)

### Experiment:

- Creation of competent cells with *B. subtilis* 168 and *B. subtilis* subsp. Natto BEST195
- Column purification and NanoDrop
- Observing results of yesterday's fibrin/milk plate assay

### Results:

- NanoDrop DNA concentration measurement for the purified PCR column was 141.8ng/μL, 97.7ng/μL, 135.8ng/μL, 133.1ng/μL for tubes 1, 2, 3, and 4 respectively
- We have determined an optimal recipe/protocol for plate creation, as well as how the strains would look depending on whether it is grown on LB medium or soybeans.





**Additional Notes:**

- Sent column samples for sequencing analysis

## 2024/08/10 (Sat)

13:00-18:00 (JST)

**Experiment Supervisor:**

Kazuyuki Fujimitsu

**Participants:**

- Shiori Kajikawa (all)

**Experiment:**

- LB agar plate creation
- NanoDrop of miniprep plasmid
- Transformation of pNK1 #4, 8, 14, 19 to E. coli C600

**Results:**

- NanoDrop DNA concentration measurement of the miniprep plasmid was 31.6ng/μL, 31.3ng/μL, 32.9ng/μL, 45ng/μL for GTstr13, 14, 15, 16 respectively.

## 2024/08/12 (Mon)

13:30-17:30 (JST)

**Experiment Supervisor:**

Fujimitsu Kazuyuki

**Participants:**

- Rikuto Fukushima (all)
- Misaki Ozawa (-16:00)
- Lee Doria (-17:00)

**Experiment:**


- Inoculation of E. coli C600
- PCR of pNK1 plasmid to create pNK2 and pNK3

	Tube		
Template DNA		1	240807_miniprep pNK1#4
		2	240807_miniprep pNK1#4

PCR program		
predenaturatio	98°C, 3 min	130 cycles
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	
extension	68°C, 30sec	
final extension	68°C, 7 min	
hold	12°C	

### Results:

- Bands were found near 1.6kbp, when it is supposed to be near 6 kbp....
- After troubleshooting, we found that the primer GT10 had another binding site, other than the one that we anticipated



### Additional Notes:

- Reconstructed primers and ordered them (GT34). Will redo this experiment on 8/21

2024/08/19 (Mon)

10:30-17:00 (JST)

### Experiment Supervisor:

Koichi Yano



**Participants:**

- Shiori Kajikawa (all)
- Mizuho Sakai (15:00-)

**Experiment:**

- Picture record of inoculation plates from 08/10
- LB plate creation (spectinomycin, chloramphenicol included)

Colonies		
	240805	C600/<05>
		C600/<110>
		C600/< $\alpha$ 1>
	240810	C600/pNK1 #4
		C600/pNK1 #8
		C600/pNK1 #14
		C600/pNK1 #19
		YAN17632

- Inoculation of E. coli C600 with various plasmids (<05>, <110>, < $\alpha$ 1>, pNK #4, 8, 14, 19)

**Results:**

None

**Additional Notes:**

## 2024/08/20 (Tue)

10:00-19:30 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Doria Lee (-12:30)
- Shiori Kajikawa (11:30-)

**Experiment:**

- Miniprep of E. coli strains that were inoculated yesterday
- NanoDrop of strains
- gDNA purification of B. subtilis (stock YAN17632)
- Creation of minimal medium for conjugational transfer
- Checking Sequence analysis of E. coli C600/pNK1 #4

**Results:**

- NanoDrop DNA concentration was 58.3ng/μL, 55.4ng/μL, 54.8ng/μL, 70.6ng/μL, 75.6ng/μL, 70.8ng/μL, 93.6ng/μL for strains GTstr17, GTstr18, GTstr19, GTstr20, GTstr21, GTstr22, and GTstr23, respectively.
- NanoDrop DNA concentration for the gDNA purification was 116.1ng/μL

**Additional Notes:**

2024/08/21 (Wed)

12:00-19:30 (JST)

**Experiment Supervisor:**

Koichi Yano


**Participants:**

- Shiori Kajikawa (all)
- Saki Tsuchiya (all)
- Mizuho Sakai (14:00-16:00)

**Experiment:**


- Electrophoresis of miniprep plasmids/B. subtilis\_YAN17632 gDNA
- Nanodrop of gDNA\_YAN17632\_240820
- PCR of pNK1 to create pNK2 and pNK3, electrophoresis
- Transformation of dimer plasmid into B. subtilis 168/pLS20cat and B. subtilis natto BEST195 (continued from 240809)

**Results:**



Each miniprep plasmid is the same length. It was divided into a trimer(9.42kbp), a dimer(5.81kbp), and a monomer(3.52kbp). The length of gDNA was 23.13kbp. Nanodrop DNA concentration measurement gDNA\_YAN17632\_240820 was 75.3ng/μL

The length of the PCR product was 5.65kbp.



#### Additional Notes:

When measuring DNA concentration with NanoDrop, we used TE buffer near NanoDrop, and it decreased a lot from yesterday. Therefore, we redid with TE buffer of our lot.

# 2024/08/22 (Thu)

11:00-17:30 (JST)

## Experiment Supervisor:

Koichi Yano

## Participants:

- Shoya Inoue (all)


## Experiment:

- Column purification of PCR samples from 08/21
- Electrophoresis of column purification sample
- Recording of plasmid transformation into *B. subtilis*/*B. subtilis* natto (photo&cell count)
- Colony PCR

Tube No.	Template
1-4	168/pL530cat /pNK1 #4(C600) colonies
5-8	168/pL530cat /<DS>(C600) colonies
9-12	168/pL530cat /noDNA colonies
13-16	168/pL530cat /KAN17632 colonies
17-20	168/pL530cat /pNK1 #4(DH5α) colonies
21-24	BEST195/pNK1 #4(C600) colonies
25-28	BEST195/<DS>(C600) colonies
29-32	BEST195/noDNA colonies
33-36	BEST195/KAN17632 colonies
37	Plasmid: pNK1 #4(C600)_240820 dil. 1/58 (1ng/ul)
38	Plasmid: <DS>(C600)_240820 dil. 1/75 (1ng/ul)


PCR program		
pre-denaturation	95°C, 5 min	
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	130 cycles
extension	68°C, 60sec	
final extension	68°C, 7 min	
hold	12°C	



## Results:



What to add	Volume plated (ml)	Colonies counted	c.f.u. (ml <sup>-1</sup> competence cells · µg <sup>-1</sup> DNA)
pNK1 #4 (CG00) miniprep240820	0.1	51	5.1E+01
	1.1	564	5.1E+03
<DS> (CG00) miniprep240820	0.1	199	2.0E+04
	1.1	904	8.0E+04
noDNA	0.1	-	-
	1.2	10	8.3E+00
YAN17632_gDNA_240820 (1/12 dil.)	0.1	3	2.0E+04
	1.1	86	8.3E+03
pNK1 #6 (DH5α) miniprep240807	0.1	0	<1E+03
	1.1	5	4.5E+03

What to add	Volume plated (mL)	Colonies count c.f.u.	(mL <sup>-1</sup> competence cells · µg <sup>-1</sup> DNA)
pNK1 #4 (CG00) miniprep240820	1.2	285	2.4E+03
<DS> (CG00) miniprep240820	1.2	330	2.8E+03
noDNA	1.2	18	1.5E+01
YAN17632_gDNA_240820 (1/32 dil.)	1.2	373	3.1E+04





**Additional Notes:**

2024/08/23 (Fri)

12:00-21:30 (JST)

### Experiment Supervisor:

Koichi Yano

### Participants:

- Shiori Kajikawa (all)

### Experiment:

- 0.8% TAE-EtBr agarose gel creation
- Electrophoresis of yesterday's colony PCR
- Quantus Fluorometer\_240822 column

The screenshot displays the QIAcuity software interface with two main sections: 'Sample Preparation' and 'DNA Quantification'.

**Sample Preparation:**

- 1/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 2/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 3/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 4/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 5/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 6/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 7/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 8/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 9/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 10/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 11/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 12/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 13/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 14/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 15/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 16/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 17/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 18/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 19/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.
- 20/206 plasmid prep:** Shows a total volume of 104.5 µl, with 110 µl of 4 hours rotation indicated in red.


**DNA Quantification:**

- 1/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 2/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 3/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 4/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 5/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 6/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 7/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 8/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 9/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 10/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 11/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 12/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 13/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 14/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 15/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 16/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 17/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 18/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 19/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.
- 20/206 plasmid prep:** Shows a concentration of 1.00 ng/µl.

- Infusion cloning to create pNK2 and pNK3







#### Additional Notes:

2024/08/24 (Sat)

11:00-16:00 (JST)

#### Experiment Supervisor:

Kazuyuki Fujimitsu

#### Participants:

- Kei Hato (11:00-14:00)
- Misaki Ozawa (12:30-)

#### Experiment:

- glycerol stock preparation
- transformation (pNK2, pNK3 to DH5α)

#### Results:

Bacterial strains culture OD value		
Sample	GTstr No.	OD(600nm)
B. subtilis 168/pLS20cat pNK1 #1	25	0.42
B. subtilis 168/pLS20cat <05> #5	26	0.64
B. subtilis 168 aprE <sup>-</sup> spc lacI/pLS20cat #	27	0.94
BEST195/pNK1 #23	28	- (error)
BEST195/<05> #26	29	1.09
BEST195 aprE <sup>-</sup> spc lacI #33	30	- (error)

#### Additional Notes:

2024/08/26 (Mon)

10:00- 19:00 (JST)

#### Experiment Supervisor:

Koichi Yano

#### Participants:

- Shiori Kajikawa (all)
- Doria Lee (-12:30)
- Misaki Ozawa (15:30-)

#### Experiment:

- LB agar plate creation
- transformation (pNK2, pNK3 to DH5α)

Tube	DNA(1μl)	DH5α
1	infusion product for pNK2 (Mix 1_240823)	25
2	infusion product for pNK3 (Mix 2_240823)	25
3	pNK1 plasmid (C600)	25
4	no DNA	25

Sp100_LB plates	Tube	Amount
240826 DH5α/pNK2-infusion 25μlまいた	1	25 ml
240826 DH5α/pNK2-infusion 残りまいた	1	500 ml
240826 DH5α/pNK3-infusion 25μlまいた	2	25 ml
240826 DH5α/pNK3-infusion 残りまいた	2	500 ml
240826 DH5α/pNK1(C600) 25μlまいた	3	25 ml
240826 DH5α no DNA 全部まいた	4	525 ml

#### Results:

- No colonies were found from the transformation from yesterday

#### Additional Notes:

• As wrong tubes were used previously for the transformation, the experiment was redone (have to use 1.5mL tubes instead of 0.6mL tubes to fit into the block) can not be heated properly when transforming

## 2024/08/27 (Tue)

10:00- 12:30 (JST)

### Experiment Supervisor:

Koichi Yano

### Participants:

• Doria Lee (all)

### Experiment:

- counting colonies from Transformation\_pNK2 into E. coli DH5α 240824
- PCR of colonies from transformation with 2 types of primer to create pNK4

Tube No.	Template
1-4	pNK3 clone 240824 TF
5-9	pNK3 clone 240826 TF
10	pNK1 colony 240826 TF
11	pNK1 (DH5α) 240807 miniprep (1ng/μl) 1 μl

PCR program		
predenaturatio	98°C, 3 min	130 cycles
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	
extension	68°C, 30sec	
final extension	68°C, 7 min	
hold	12°C	

Tube No.	Template
21-22	pNK3 clone 240824 TF #1-2
23-24	pNK3 clone 240826 TF #5-6

PCR program		
predenaturatio	98°C, 3 min	130 cycles
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	
extension	68°C, 1min10sec	
final extension	68°C, 7 min	
hold	12°C	

## Results:

Transformation_pNK2 into E. coli DH5α 240824 (Plate:SP 100)		
Plates 1-6		Colonies counted
1	240824 DH5α/pNK2-infusion 25μlまいた	0
2	240824 DH5α/pNK2-infusion 残りまいた	0
3	240824 DH5α/pNK3-infusion 25μlまいた	0
4	240824 DH5α/pNK3-infusion 残りまいた	4
5	240824 DH5α/pNK1(C600) 25μlまいた	60
6	240824 DH5α no DNA 全部まいた	0

Transformation_pNK2 into E. coli DH5α 240826 *USING 1.5mL Tubes (Plate:SP 100)		
Plates 1-6		Colonies counted
1	240826 DH5α/pNK2-infusion 25μlまいた	0
2	240826 DH5α/pNK2-infusion 残りまいた	0
3	240826 DH5α/pNK3-infusion 25μlまいた	0
4	240826 DH5α/pNK3-infusion 残りまいた	5
5	240826 DH5α/pNK1(C600) 25μlまいた	1
6	240826 DH5α no DNA 全部まいた	0

## Additional Notes:

2024/08/28 (Wed)

10:30- 16:00 (JST)

## Experiment Supervisor:

Kazuyuki Fujimitsu

## Participants:

- Kei Hato (-14:00)
- Misaki Ozawa (11:00 - )



## Experiment:

- Electrophoresis of yesterday's colony PCR

Tube No.	Template
1-4	pNK3 clone 240824 TF
5-9	pNK3 clone 240826 TF
10	pNK1 colony 240826 TF
11	pNK1 (DH5α) 240807 miniprep (1ng/μl) 1 μL

- Inoculation of pNK3 from electrophoresis
- Planting colonies (BEST 195 and B. subtilis 168) to difco and star agar plates

Results:



Additional Notes:

## 2024/09/02 (Mon)

10:00- 13:00 (JST)

### Experiment Supervisor:

Kazuyuki Fujimitsu

### Participants:

• Kei Hato (all)

### Experiment:

- Glycerol freeze stock of bacterial strain DH5a/pNK3 clone, miniprep
- Column purification of PCR sample from 8/27

DNA sample	Sample data:
240827_PCR_#21	240828
240827_PCR_#24	240828

### Results:

None

### Additional Notes:

## 2024/09/03 (Tue)

16:30-19:00 (JST)

### Experiment Supervisor:

Koichi Yano

### Participants:

• Wingdor Doria Lee (all)

### Experiment:

- Quantus of Fluoromenter\_240822 column
- In-Fusion of Fluoromenter\_240822 column
- Miniprep of E. coli DH5α/pNK3 #6, inoculation into LB+Sp for glycerol stock creation

### Results:

None

### Additional Notes:

# 2024/09/04 (Wed)

16:30-19:00 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Wingdor Doria Lee (all)

**Experiment:**

- LB agar plate creation
- Glycerol Freeze Stock of bacterial strains

**Results:**

None

**Additional Notes:**

# 2024/09/05 (Thu)

14:00-19:30 (JST)

**Experiment Supervisor:**

Koichi Yano

**Participants:**

- Mizuho Sakai (all)
- Wingdor Doria Lee (16:00-19:00)

**Experiment:**

- Transformation\_pNK5 into E. coli DH5α 240904 \*USING 1.5mL Tubes

3. Defreeze E. coli DH5α on ice and mix thoroughly

1. Mix 15ul of E. coli DH5α competent cell with 5ul of plasmid DNA in 1.5 ml tube, labeled as below

Tube	Volume (ul)	Label
3. Infusion product for pNK5 (DNA 1, 240904)	25	
3. no DNA	25	

2. Put on ice for 10 min, incubate at 42°C for 45 sec, put on ice for 2 min

3. Add 500 ul of LB medium, put in 17 degrees for 30 min. Aliquot 25ul of tube 1, 2 into LB plate.

Before aliquot, take 200ul of fresh cell residues onto the plate, so that it won't stop in too fast.

4. Count the remaining tube foods, 10,000 spm for 2 min. Freeze 100ul of supernatant, then wash the rest, and resuspend cells with the preserved sup.

5. Spread the remaining liquid onto LB plate

For steps 3 and 5, reference the chart below


Sp100 LB plates	Tube	Amount
242904 (E. coli DH5α) medium (17.4°C)	1	25 ul
240904 (E. coli DH5α) medium (17.4°C)	1	500 ul
242904 (E. coli DH5α) medium (17.4°C)	2	525 ul

6. incubate at 37 degrees, O/N



- QIAprep Spin Miniprep Kit (pNK3 #6)
- Electrophoresis (pNK3 #1, 6, 9, 21, 24)

### Results:



### Additional Notes:

2024/09/08 (Sun)

11:00-18:30 (JST)

### Experiment Supervisor:

Kazuyuki Fujimittsu

### Participants:

- Misaki Ozawa (all)


### Experiment:

- PCR and electrophoresis of colonies from DH5a/pNK5-infusion

Tube No.	Template	
1-3	DH5a/pNK5-infusion 240904 #1	
4-6	DH5a/pNK5-infusion 240904 #2	
7-9	DH5a/pNK5-infusion 240904 #3	
10-12	DH5a/pNK5-infusion 240904 #4	
PCR program		
predenaturation	98°C, 3 min	730 cycles
denaturation	98°C, 10 sec	
annealing	55°C, 5 sec	
extension	68°C, 1min10sec	
final extension	68°C, 7 min	
hold	12°C	

## Results:

Lane No.	Tube #	10X	Ex Loading Dye	TE	Ladder (PCR)	Total	
1		-	-	-	-	5	µL
2	1	1	1	4	-	5	µL
3	2	1	1	4	-	5	µL
4	3	1	1	4	-	5	µL
5	4	1	1	4	-	5	µL
6	5	1	1	4	-	5	µL
7	6	1	1	4	-	5	µL
8	7	1	1	4	-	5	µL
9	8	1	1	4	-	5	µL
10	9	1	1	4	-	5	µL
11	10	1	1	4	-	5	µL
12	11	1	1	4	-	5	µL
13	-	-	-	5	5	5	µL
14	12	1	1	4	-	5	µL



## Additional Notes:

- Used the wrong ladder, so redid electrophoresis a second time.

# 2024/09/10 (Mon)

16:00-18:00 (JST)

## Experiment Supervisor:

Koichi Yano

## Participants:

• Shoya Inoue (all)

## Experiment:

• Combine DNA samples from 9/8 colony PCR

	DNA sample	Sample data:
A	240908_colonyPCR #1~3 mix	240908
B	240908_colonyPCR #4~6 mix	240908
C	240908_colonyPCR #7~9 mix	240908
D	240908_colonyPCR #10~12 mix	240908

## Results:

None

## Additional Notes:

# 2024/09/11 (Wed)

17:00-20:00 (JST)

## Experiment Supervisor:

Koichi Yano

## Participants:

• Misaki Ozawa (all)

## Experiment:

• Glycerol Freeze Stock of DH5a/pNK5-infusion  
• Miniprep of DH5a/pNK5-infusion

## Results:

None

## Additional Notes:

# 2024/09/16 (Mon)

11:00-18:30 (JST)

## Experiment Supervisor:

Kazuyuki Fujimitsu

## Participants:

- Shoya Inoue (all)
- Mizuho Sakai (13:00-)

## Experiment:

- Transformation\_pNK2 into E. coli DH5α 240826
- NanoDrop DNA conc. measurement(240910\_column)

## Results:

- NanoDrop DNA concentration values for the four samples from 9/10 was 27.7ng/μL, 34.5ng/μL, 90.3ng/μL, and 79.5ng/μL
- When spreading on LB medium, mistakenly used an item that contained Cm.

## Additional Notes:

# 2024/09/17 (Tue)

12:00-19:30 (JST)

## Experiment Supervisor:

Koichi Yano

## Participants:

- Saki Tsuchiya (all)
- Mizuho Sakai (16:30-)

## Experiment:

- Sequence analysis

240917		シーケンス実験_02EM		サンプル測定 遺伝子					
No.	Date	Sample No.	DNA	Primer	DNA amount (ng)		24h濃度 H2O (μL)	備考	
					Primer 1 (μL) 1.00	DNA conc (ng/μL)			
1	240917	1	240910_column_A	2716	1	27.7	1.01	96.39	240917_#1
2	240917	2	240910_column_B	2716	1	34.5	1.39	17.70	240917_#2
3	240917	3	240910_column_C	2716	1	90.3	1.11	98.88	240917_#3
4	240917	4	240910_column_D	2716	1	79.5	1.28	98.74	240917_#4

- Electrophoresis (100V, 25min, 2%agarose gel )
- Transformation\_pNK2 into E. coli DH5α 240826

Tube	DNA(μl)		C600
1 C600-1	pNK5 #1_240911		100
2 C600-2	pNK5 #2_240911		100
3 C600-3	pNK5 #3_240911		100
4 C600-4	pNK5 #4_240911		100
5 C600-ND	-		100

つまりこのグリセロールには何もいれない。ネガコン。

Sp100 LB plates	Tube	Amount
240917 C600/pNK5 #1_240911 全部まいた	1	525 μl
240917 C600/pNK5 #2_240911 全部まいた	2	525 μl
240917 C600/pNK5 #3_240911 全部まいた	3	525 μl
240917 C600/pNK5 #4_240911 全部まいた	4	525 μl
240917 C600 no DNA 全部まいた	5	525 μl


- Transformation of dimer plasmid into B. subtilis 168/pLS20cat (continued from 240809)

Tube	What to add	How much (μl)
1	pNK5 (DH5α) #1_240911	20
2	pNK5 (DH5α) #2_240911	20
3	pNK5 (DH5α) #3_240911	20
4	pNK5 (DH5α) #4_240911	20
5	no DNA	0

- Creating Natto with soybean and NK Day1

## Results:

- In electrophoresis, all bands were seen. Only pnc5 3 was seen in different lengths.



## Additional Notes:

We redid the "Transformation\_pNK2 into E. coli DH5α 240826" process that I made a mistake in yesterday.

# 2024/09/18 (Wed)

17:00-21:00 (JST)

## Experiment Supervisor:

•Koichi Yano

## Participants:

•Arisa Tani (all)

## Experiment:

- Creating Natto with soybean and NK Day2
- Colony PCR (168/pLS20cat pNK5)

Tube No.	Template		
1 - 4	168/pLS20cat /pNK5 #1 colonies		
5 - 9	168/pLS20cat /pNK5 #2 colonies		
10 - 14	168/pLS20cat /pNK5 #3 colonies		
15 - 19	168/pLS20cat /pNK5 #4 colonies		
20	Positive control (plasmid pNK5 #1)		
21	Positive control (plasmid pNK5 #3)		
PCR program			
predenaturation	98°C	3 mins	x 30
denaturation	98°C	10 secs	
annealing	55°C	5 secs	
extension	68°C	30 secs	
final extension	68°C	7 mins	
hold	12°C	-	

- Culture of transformed 168, (not transformed) BEST195 and S903 in preparation for conjugal Transmission

Donor	Str No.	Antibiotics	Recipient	Str No.	Antibiotics
168/pLS20cat pNK1	25	Sp100 Cm5	BEST195	1	None
168/pLS20cat <05>	26	Sp100 Cm5	S903	35	None
168/pLS20cat pNK5 #1 - 1		Sp100 Cm5			
168/pLS20cat pNK5 #1 - 2		Sp100 Cm5			
168/pLS20cat pNK5 #2 - 5		Sp100 Cm5			
168/pLS20cat pNK5 #2 - 6		Sp100 Cm5			
168/pLS20cat pNK5 #3 - 10		Sp100 Cm5			
168/pLS20cat pNK5 #3 - 11		Sp100 Cm5			
168/pLS20cat pNK5 #4 - 15		Sp100 Cm5			
168/pLS20cat pNK5 #4 - 16		Sp100 Cm5			

## Results:

None

## Additional Notes:

# 2024/09/19 (Thu)

16:00-20:30 (JST)

**Experiment Supervisor:**

- Koichi Yano

**Participants:**

- Rikuto Fukushima (-18:00)

- Arisa Tani (17:00-)

**Experiment:**

- Milk plate assay

- Fibrin plate assay

- Creating Natto with soybean and NK Day3

- Conjugation DAY 1

- Preparation of Glycerol Freeze Stock

**Results:**

Natto creation (sample) is completed.



## 2024/09/20 (Fri)

17:00-19:00 (JST)

**Experiment Supervisor:**

- Koichi Yano

**Participants:**

- Misaki Ozawa (all)

**Experiment:**

- Glycerol Freeze Stock of bacterial strains
- Conjugation Day2

**Results:**

None

**Additional Notes:**

2024/09/22 (Sun)

11:00-17:30 (JST)

**Experiment Supervisor:**

Kazuyuki Fujimitsu

**Participants:**

- Arisa Tani (all)

**Experiment:**

- S903, BEST195 fibrin plate assay preliminary experiment
- Colony PCR


Number	Plasmid / Strain	Samples	Primer_F	Primer_R
1	pNK1 / 195	4	GT16	GT18
2	<05> / 195	4	GT16	GT17
3	pNK5 (O/N) / 195 (O/N)	4	GT16	GT18
4	pNK5 / 195	4	GT16	GT18
5	pNK1 / S903	4	GT16	GT18
6	<05> / S903	4	GT16	GT17
7	pNK5 / S903	4	GT16	GT18
8	Positive control (pNK5-#1)	1	GT16	GT18

PCR program			
predenaturation	98°C	3 mins	x 30
denaturation	98°C	10 secs	
annealing	55°C	5 secs	
extension	68°C	30 secs	
final extension	68°C	7 mins	
hold	12°C	-	

- Electrophoresis of colony PCR product

**Results:**





## Yesterday's Conjugation

Results of conjugation						
MM + Ca	Donor / Recipient	BEST195	5903	MM + Ca (overnight)	Donor / Recipient	BEST195
	168/pLS20cat pNK1	>2,000	>2,000		168/pLS20cat pNK5 #1 - 2	16
	168/pLS20cat <05>	>2,000	>2,000			
	168/pLS20cat pNK5 #1 - 2		6			35
MM - Ca	Donor / Recipient	BEST195	5903	MM - Ca (overnight)	Donor / Recipient	BEST195
	168/pLS20cat pNK1	>2,000	>2,000		168/pLS20cat pNK5 #1 - 2	5
	168/pLS20cat <05>	>2,000	>2,000			
	168/pLS20cat pNK5 #1 - 2		1			19

• The colony is clearly visible only on the plate that transmitted NK5. This is unknown because we did not take a negative control, but it may appear spontaneously. If it is spontaneous, the transmission of NK5 failed, while the other >2,000 colonies, pNK1 and <05>, were successfully transmitted. The successful plate has the potential to have produced >2,000 colonies because too many were spread on the plate.

## Additional Notes:

2024/09/23 (Mon)

13:00-20:00 (JST)

## Experiment Supervisor:

Kazuyuki Fujimitsu

## Participants:

- Arisa Tani (all)

### Experiment:

- Overnight culture creation
- Culture creation for PCR
- LB agar plates for filter membrane method
- Minimal medium agar plates

### Results:

None

### Additional Notes:

## 2024/09/24 (Tue)

13:30-22:00 (JST)

### Experiment Supervisor:

Koichi Yano

Kazuyuki Fujimitsu

### Participants:

- Mizuho Sakai (all)
- Saki Tsuchiya (17:00-20:00)

### Experiment:

- Preparation of LB medium for Colony PCR
- Colony PCR

Number	Plasmid / Strain	Colony	Samples	Primer_F	Primer_R
4	pNK1 / 195	1-1	4	GT18	GT18
2	pNK1 / 195	1-3	4	GT16	GT18
3	<O5> / 195	2-3	4	GT16	GT17
4	pNK5 (G/N) / 195 (G/N)	3-4	4	GT16	GT18
5	pNK5 / S003	7-3	4	GT16	GT18
6	pNK5 / S003	7-4	4	GT16	GT17
7	Positive control (pNK5-#1)		1	GT16	GT18
8	Positive control (<O5>)		1	GT16	GT17

predenaturation	98°C	3 mins	x 30
denaturation	98°C	10 secs	
annealing	55°C	5 secs	
extension	68°C	30 secs	
final extension	68°C	7 mins	
hold	12°C	-	

- Confirmation of the bacteria (Bacillus subtilis or Bacillus subtilis natto)
- Preparation of NB medium for fibrin assay

**Results:**

None

**Additional Notes:**

## 2024/09/25 (Wed)

17:00-22:00 (JST)

**Experiment Supervisor:**

Koichi Yano

Kazuyuki Fujimitsu

**Participants:**

- Arisa Tani (all)

**Experiment:**

- Conjugation Day2

- Inoculated on LB medium

- Electrophoresis (yesterday's Colony PCR)

**Results:**

OD value measured before inoculating LB medium

Strain	OD	amount for 10 <sup>9</sup>
BEST195	0.91	1.373626374
S903	1.35	0.9259259259
pNK1	1.7	0.7352941176
<05>	0.95	1.315789474
N.C.	1.41	0.8865248227
pNK5	2.19	0.5707762557

- Electrophoresis (yesterday's Colony PCR)



#### Additional Notes:

2024/09/26 (Thu)

17:00-22:00 (JST)

#### Experiment Supervisor:

Koichi Yano

Kazuyuki Fujimitsu

#### Participants:

- Arisa Tani (all)
- Mizuho Sakai (all)

#### Experiment:

- Preparation of 100mL of NB medium, fibrin assay

Nutrient Broth		
powder	0.8g	
DW	100mL	x 2 bottles

- Colony PCR, Colony count


Number	Plasmid / Strain	Colony No.	Primer F	Primer R
1	BEST285/-OG-	1		GT17
2		2		
3		3		
4	BEST285/pM4	4		GT18
5		5		
6		6		
7		7		
8		8		
9		9		
10	BEST285/pM4	10		GT18
11		11		
12		12		
13	BEST285/pM4	13		GT17
14		14		
15		15		
16	S903/OG-	16	GT18	GT17
17		17		
18		18		
19		19		
20	S903/pM4	20		GT18
21		21		
22		22		
23		23		
24	S903/pM4	24		GT18
25		25		
26	S903/pM4	26		GT17
27	S903/pM4	27		GT17
28	S903/pM4	28		GT18
29	S903/pM4	29		GT18

PCR Cycle			
predenaturation	98°C	3 mins	
denaturation	98°C	10 secs	x 30
annealing	55°C	5 secs	
extension	68°C	30 secs	
final extension	68°C	7 mins	
hold	12°C	-	

- Electrophoresis of Colony PCR products
- fibrin plate assay

## Results:

- Electrophoresis of Colony PCR products




Colony count

50µl			50µl		
Donor / Recipient	BEST195	S903	Donor / Recipient	BEST195	S903
168/pL520cat <05>	7	0	168/pL520cat <05>	6	7
168/pL520cat pNR1	1	0	168/pL520cat pNR1	30	8
168/pL520cat pNR5 #3-15	11	0	168/pL520cat pNR5 #3-15	41	48
168/pL520cat (N.C.)	>3,000*	1	168/pL520cat (N.C.)	62	>3,000*


\*The colonies' sizes are too small to count.


## Conjugation results

### BEST195



### S903





**Additional Notes:**

2024/09/27 (Fri)

12:30-16:30, 17:30-22:00 (JST)

**Experiment Supervisor:**

Koichi Yano

Kazuyuki Fujimitsu

**Participants:**

• Arisa Tani (17:30-22:00)


• Mizuho Sakai (12:30-16:30)

**Experiment:**

- Confirmation of results of fibrin assay
- Preparation of NB medium
- Preparation of Conjugation

**Results:**

Results of yesterday's fibrin plate assay



#### **Additional Notes:**

- The bacteria did not increase in the NB medium made in yesterday's "Preparation of NB medium" process. The use of thin tubes may have resulted in a lack of oxygen due to the small surface area.

2024/09/28 (Sat)

13:00-21:30 (JST)

#### **Experiment Supervisor:**

Kazuyuki Fujimitsu

#### **Participants:**

- Mizuho Sakai (all)
- Arisa Tani (14:00-)

#### **Experiment:**

- Conjugation
- fibrin plates creation
- Preparation NB medium for fibrin assay, fibrin assay
- Fibrin plate assay
- Made some glycerol stocks

#### **Results:**

OD value of inoculated NB medium for fibrin assay



Number	plasmid/strain	OD	dillution for OD=0.5
3	BEST195 (GTstr.1)	1.17	0.234
4	BEST195/pNK1 (3)	1.17	0.234
5	BEST195/pNK5 (3)	2.16	0.432
6	S903/<05> (3)	0.94	0.188
7	S903/pNK1 (4)	1.53	0.306
8	S903/pNK5 (1)	1.83	0.366
9	BEST195/<05>	1.41	0.282
10	S903 (GTstr.35)	1.2	0.24

• Made some glycerol stocks

Glycerol stock	
4 BEST195/pNK1 (3)	47
5 BEST195/pNK5 (3)	48
9 BEST195/<05>	49
6 S903/<05> (3)	50
7 S903/pNK1 (4)	51
8 S903/pNK5 (1)	52

**Additional Notes:**

2024/09/29 (Sun)

11:00-22:00 (JST)

**Experiment Supervisor:**

Kazuyuki Fujimitsu

**Participants:**

• Arisa Tani (all)

**Experiment:**

• Creation of sodium phosphate buffer

Sodium phosphate buffer		
NaH <sub>2</sub> PO <sub>4</sub> 2H <sub>2</sub> O	1.919	g
Na <sub>2</sub> HPO <sub>4</sub> 12H <sub>2</sub> O	13.503	g
Water	400+x	mL
Total	500	mL

• Fibrin plate creation


• Halo assay with natto suspension liquid and LB medium

**Results:**


• Yesterday's results of fibrin plate assay

16h


BEST195×3



BEST195×4




S903×3




20h


BEST195×3




BEST195×4




S903×3




24h  
BEST195×3



BEST195x4



S903x3



• OD value of natto solution in the process of “Halo assay with natto suspension liquid.

Strain	OD	dillution for OD=0.5
1 BEST195	2.92	0.0584
2 BEST195/pNK1	4.64	0.0928
3 BEST195/pNK5	6.08	0.1216
4 BEST195/<05>	6	0.12
5 S903	2.78	0.0556
6 S903/pNK1	4.16	0.0832
7 S903/pNK5	3.92	0.0784
8 S903/<05>	1.74	0.0348

#### Additional Notes:

## 2024/09/30 (Mon)

18:00-22:00 (JST)

#### Experiment Supervisor:

Koichi Yano

Kazuyuki Fujimitsu

#### Participants:

- Mizuho Sakai (all)
- Arisa Tani (19:00-)

#### Experiment:

- Fibrin plate assay
  - Liquid inoculated on LB medium
  - Liquid suspended natto

#### Results:

None

## 2024/10/01 (Tue)

13:30-22:00 (JST)

#### Experiment Supervisor:

Koichi Yano

Kazuyuki Fujimitsu

#### Participants:

- Mizuho Sakai (all)



• Arisa Tani (19:00-)

## Experiment:

- Confirmation of yesterday's results of the fibrin assay
- Colony PCR

Number	Plasmid / Strain	Primer_F	Primer_R	Amplified length (bp)
1	PC-198/ΔOS (LB culture (centrifuge culture, pellet direct into PCR))	GT16	GT17	1,588
2	ΔOS>RC- (240820 G10r21 C600/ΔOS) 1/75 dil. 3μL			1,588
3	BEST195/ΔOS (240927 NB culture No.3, centrifuge pellet, suspend directly into PCR mix)			1,588
4	BEST195/ΔOS (240927 NB culture No.3, centrifuge pellet, suspend directly into PCR mix)			1,588

PCR Cycle				
predenaturation	90°C	3 mins		
denaturation	98°C	10 secs		
annealing	55°C	5 secs	x 30	
extension	68°C	30 secs		
final extension	68°C	7 mins		
hold	12°C	-		

- Electrophoresis of Colony PCR products


## Results:

- Yesterday's results of fibrin assay


**Liquid inoculated on LB medium**

**16h**


**BEST195**




**S903**




20h  
BEST195




S903



24h  
BEST195




S903




Natto suspension liquid

16h

BEST195




S903




20h

BEST195







S903



24h  
BEST195




S903



- OD values measured before Colony PCR

	OD
BEST195/05 (240927 NB culture No.2, centrifuge pellet, suspend directly into PCR rxn)	0.92
BEST195/05 (240927 NB culture No.9, centrifuge pellet, suspend directly into PCR rxn)	1.31

- Electrophoresis



#### Additional Notes:

There was a little water on the surface of the fibrin plate. We were able to confirm the halo, but need to keep a close eye on it.